

News

Clinical Laboratory Test Update

Coagulation Testing Reference Range Updates

Effective September 4, 2018, the University of Colorado Hospital Clinical Laboratory has revised the reference ranges for the following tests to reflect the sensitivity of the test methods.

TEST	EPIC CODE	Old Reference Range	New Reference Range
Prothrombin Time	LAB320	11.8 – 14.3 seconds	11.8 – 13.8 seconds
INR	LAB320	0.9 – 1.1	0.9 – 1.1 (no change)
Activated Partial Thromboplastin Time (APTT)	LAB325	23.9 – 34.1 seconds	23.8 – 34.4 seconds
Therapeutic Range of APTT for Heparin Monitoring		80.7 – 122.7 seconds	80.7 – 122.7 seconds (no change)
Activated Partial Thromboplastin Time, Heparin Neutralized (APTT, Heparin Neutralized) - <i>only used in rare occasion when sample cannot be recollected and there is suspected contamination of heparin</i>	LAB4315	23.9 – 34.1 seconds	23.8 – 34.4 seconds (no change)

Please call Stuart Lind, MD at 303-724-1383 if you have any questions or visit our website at <https://www.uchealth.org/professionals/Pages/Clinical-Laboratory.aspx> for additional information.

Issued: August 15, 2018