

News

Clinical Laboratory Test Update

SARS CoV-2 IgG (nucleocapsid) antibody testing update

Effective June 13, 2022 the University of Colorado Hospital Clinical Laboratory will discontinue SARS CoV-2 IgG (nucleocapsid) antibody testing. Samples received in the lab will be sent to our reference lab for testing.

Test name	COVID-19 Antibody IgG-SARS CoV-2
EAP	LAB1939
Collect	Serum separator tube (SST) or EDTA plasma
Storage/Transport	Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 0.5mL serum or plasma to an ARUP Standard Transport Tube. Minimum volume: 0.25 mL serum or plasma.
Stability	Refrigerated: 1 week; Frozen: 1 month
Method	Qualitative Chemiluminescent Immunoassay
Performed	Sunday - Saturday
Reported	1 – 5 days
Reference Interval	Negative
CPT	86769

Please call Brian Harry, MD PhD at 303-848-7070 if you have any questions or visit our test catalog at <https://www.testmenu.com/universityhospital> for additional information.

Issued: 06/15/2022