

# News

## Clinical Laboratory Test Update

### Hexagonal Phospholipid Neutralization

Effective Thursday, February 22, 2017 the reference range for Hexagonal Phospholipid Neutralization will be updated. This change will better reflect our normal patient population with the current method.

Test Mnemonic:       HEXPPOS  
Epic Code:            LAB3159  
Methodology:        Clotting assay  
CPT Code:            85598

	<del>OLD</del> Reference Range	<del>NEW</del> Reference Range
HEXPPOS	Negative <8.6 seconds	Negative <9.0 seconds

Please call Stuart Lind, MD at 303-724-1383 if you have any questions or visit our website at <https://www.uchealth.org/professionals/uch-clinical-laboratory/> for additional information.

Issued: 02/15/2017