

News

Clinical Laboratory Test Update

Beta-2-microglobulin (B2M)

Effective August 24, 2020 the University of Colorado Hospital Clinical Laboratory is sending Beta-2-microglobulin samples to Medical Center of the Rockies for testing.

Test Name	Beta-2-microglobulin serum (B2M)
Methodology	Quantitative Immunoturbidimetric
EPIC code	LAB49
Collect	One 5 mL red top tube Minimum volume: 0.5 mL serum
Stability	Internal: Deliver to lab immediately at ambient temperature Offsite: Centrifuge within 2 hours of collection and deliver to lab at 2 – 8 degrees C
Performed	Batch Testing: Monday, Wednesday, Friday
Reported	Same day of testing
CPT	82232
Reference Range	0.97 – 1.84 mg/dL

Please call Mary Ellen Boyle, at 720-848-7067 if you have any questions or visit our test directory at <https://www.testmenu.com/universityhospital> for additional information.