

News

Clinical Laboratory Test Update

TSH Reflex

Effective February 13, 2023, the TSH Reflex test will be available at Cherry Creek Medical Center Clinical Laboratory. The standard non-reflexed version of this test (TSH Ultrasensitive; LAB129) is still also available.

Please call Cherry Creek Medical Center Laboratory at 720.516.9240 if you have any questions or visit our test catalog at <https://www.testmenu.com/universityhospital> for additional information.

Test Name and Epic Order Code	TSH Reflex (LAB3787)
Collect	4.5 mL Gold or Red (SST) Also acceptable: 6 mL Plain Red, no gel
Storage/Transport	Internal: Deliver to lab immediately at ambient temperature. Offsite: Centrifuge within 30 minutes of collection and transport to lab at 2-8 degrees C.
Stability	Ambient, 8 hours; Refrigerated, 7 days; Frozen, 4 weeks.
Reflex	Abnormal TSH values outside the reference range will reflex to perform a Free T4 Test (FT4).

Issued: 2/13/2023